



X-TPLO Catalogue

X-TPLO Plates

The X-TPLO plates are available in four sizes, divided into two groups based on size: medium and large. The medium plates are generally used for dogs weighing up to 35 kg, while the large plates are intended for giant breed dogs weighing over 35 kg. All sizes are compatible with Synthes screws.



Medium



Large

Size	<u>Right</u>	<u>Left</u>
<u>Medium</u>	100020201B	100020202B
<u>Large</u>	100020203B	100020204B

Intraoperative X-TPLO guides

The intraoperative X-TPLO guides are available in four versions, each specifically designed to match a corresponding plate size and orientation: medium right, medium left, large right and large left.



Medium Right



Medium Left



Large Right



Large Left

Size	<u>Right</u>	<u>Left</u>
<u>Medium</u>	S00040017A	S00040018A
<u>Large</u>	S00040019A	S00040020A

Ratchet compression forceps



Code

S00040021A

Progressive reduction forceps



Code

S00040022A

Synthes-compatible locking screws



<i>Size</i>	<i>Code</i>
Ø 3.5 x 24 mm	V00070001A
Ø 3.5 x 26 mm	V00070002A
Ø 3.5 x 28 mm	V00070003A
Ø 3.5 x 30 mm	V00070004A
Ø 3.5 x 32 mm	V00070005A
Ø 3.5 x 34 mm	V00070006A
Ø 3.5 x 36 mm	V00070007A
Ø 3.5 x 38 mm	V00070008A
Ø 3.5 x 40 mm	V00070009A
Ø 3.5 x 45 mm	V00070010A

Synthes-compatible cortical screws



<i>Size</i>	<i>Code</i>
Ø 3.5 x 20 mm	V00070101A
Ø 3.5 x 22 mm	V00070102A
Ø 3.5 x 24 mm	V00070103A
Ø 3.5 x 26 mm	V00070104A
Ø 3.5 x 28 mm	V00070105A
Ø 3.5 x 30 mm	V00070106A

Price List

Code	Description	Price
I00020201B	Medium right X-TPLO plate, 316L	85,00 €
I00020202B	Medium left X-TPLO plate, 316L	85,00 €
I00020203B	Large right X-TPLO plate, 316L	105,00 €
I00020204B	Large left X-TPLO plate, 316L	105,00 €
S00040017A	Medium right intraoperative XTPLO guide	30,00 €
S00040018A	Medium left intraoperative XTPLO guide	30,00 €
S00040019A	Large right intraoperative XTPLO guide	30,00 €
S00040020A	Large left intraoperative XTPLO guide	30,00 €
S00040021A	Ratchet compression forceps	60,00 €
S00040022A	Progressive reduction forceps	60,00 €
V00070001A	Synthes-compatible locking screw Ø 3.5 x 24 mm	12,50 €
V00070002A	Synthes-compatible locking screw Ø 3.5 x 26 mm	12,50 €
V00070003A	Synthes-compatible locking screw Ø 3.5 x 28 mm	12,50 €
V00070004A	Synthes-compatible locking screw Ø 3.5 x 30 mm	12,50 €
V00070005A	Synthes-compatible locking screw Ø 3.5 x 32 mm	12,50 €
V00070006A	Synthes-compatible locking screw Ø 3.5 x 34 mm	12,50 €
V00070007A	Synthes-compatible locking screw Ø 3.5 x 36 mm	12,50 €
V00070008A	Synthes-compatible locking screw Ø 3.5 x 38 mm	12,50 €
V00070009A	Synthes-compatible locking screw Ø 3.5 x 40 mm	12,50 €
V00070010A	Synthes-compatible locking screw Ø 3.5 x 45 mm	12,50 €
V00070101A	Synthes-compatible cortical screw Ø 3.5 x 20 mm	11,50 €
V00070102A	Synthes-compatible cortical screw Ø 3.5 x 22 mm	11,50 €
V00070103A	Synthes-compatible cortical screw Ø 3.5 x 24 mm	11,50 €
V00070104A	Synthes-compatible cortical screw Ø 3.5 x 26 mm	11,50 €
V00070105A	Synthes-compatible cortical screw Ø 3.5 x 28 mm	11,50 €
V00070106A	Synthes-compatible cortical screw Ø 3.5 x 30 mm	11,50 €

Ad Maiora has the goal of developing safe and effective surgical techniques for veterinary patients and to create tools that facilitate the task of the surgeon. Thanks to the intellectual input and direct experience of the surgeons who cooperate for this development, Ad Maiora is constantly providing solutions in an attempt to give effective answers to the problems encountered in the surgical activity.



www.ad-maiora.eu

